Tube


Fig. 1

ST-20/50 ...

DT-20/50 ...

BMT-20/50 ...

Europe Middle East Africa

IKA® - Werke
GmbH & Co.KG
Janke & Kunkel-Str. 10
D-79219 Staufen
Tel.: +49 7633 831-0
Fax: +49 7633 831-98
E-Mail: sales@ika.de

North America

IKA® Works, Inc.
2635 North Chase
Pkwy SE
Wilmington
NC 28405-7419 USA
Tel.: 800 733-3037
Tel.: +1 910 452-7059
Fax: +1 910 452-7693
E-Mail: usa@ika.net

China

IKA® Works Guangzhou
173-175 Friendship Road
Guangzhou
Economic and Technological
Development District
510730 Guangzhou, China
Tel.: +86 20 8222-6771
Fax: +86 20 8222-6776
E-Mail: sales@ikagz.com.cn

Asia Australia

IKA® Works (Asia)
Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara
Technology Park
47810 Petaling Jaya
Selangor, Malaysia
Tel.: +60 3 7804-3322
Fax: +60 3 7804-8940
E-Mail: sales@ika.com.my

Japan

IKA® Japan K.K.
293-1 Kobayashi-cho
Yamato Konyama Shi, Nara
639-1026 Japan
Tel.: +81 743 58-4611
Fax: +81 743 58-4612
E-Mail: info@ika.ne.jp

Korea

IKA® Korea Co LTD
1710 Anyang Trade Center
1107 Buhung-dong,
Dongan-gu
Anyang City, Kyeonggi-do
Post code: 431-817
South Korea
Tel.: +82 31-380-6877
Fax: +82 31-380-6878
E-Mail: michael@ikakorea.co.kr

India

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Safety instructions EN

- Read the operating instructions in full before starting up and follow the safety instructions.**
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Keep the operating instructions in a place where they can be accessed by everyone.

- Ensure that only trained staff work with the appliance.
- Ensure that the cover is screwed on tightly.
- The tube must only be attached and removed whilst the motor is stationary.
- Ensure that the tube is firmly attached to the bayonet lock connector on the drive unit prior to operating the device.
- Ensure that the tube is able to withstand the chemical properties of the selected medium by filling the tube and storing for a suitable period time.
- Always open the tube carefully after use as the media in the tube may heat up due to transfer of energy during operation, leading to pressurization of the container: risk of material spraying, protective equipment must be worn.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles.
- The temperature of the material must not exceed 40 °C.
- Note the operating instructions of the disperser unit.

Product information

ST-20/50... comprise mixing vessel and stirring tool, DT-20/50... comprise mixing vessel and dispersion tool, BMT-20/50... comprise mixing vessel, stirring tool and balls.

ST-20, ST 20-M, ST-50, ST-50-M, DT-20, DT 20-M, DT-50, DT-50-M and BMT-20, BMT 20-M, BMT-50, BMT-50-M are not sterile and not free of DNA or RNA.

ST-20-gamma, ST-20-M-gamma,DT-20-gamma, DT-20-M-gamma and BMT-20-gamma, BMT-20-M-gamma are sterilized.

TC-20-M und TC-50-M are tube covers with a pierceable membran for needles or pipette points.

The tubes are ideal for applications, with which no cross contamination is permissible.

The tubes ST-20, DT-20 und BMT-20 are appropriate only for a single-use. Multiple-use can destroy the tube and the disperser unit.

Application instruction

Store the tubes in a dry place, in their packaging and at room temperature.

Note the use-by date of the tubes.

Protect the tubes against ultraviolet light.

Attention! Never run the DT-20/50... dry.

For ST 20...,DT-20... and BMT-20... the fill level of the medium must be between the 2 respectively 5 ml and max. 15 ml marks. For ST 50..., DT-50... und BMT-50... the fill level of the medium must be between the 15 ml and max. 50 ml marks. Don't turn the tube after piercing the membran, then the leak-tightness is ensured no longer.

Always position the sample in the centre of the dispersing unit.

The solid bodies in samples must not be frozen or too hard.

Dispersing causes the medium to heat up.

Likewise the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute.

Longer application times bring no improvements, increase however the sample temperature substantially.

Dispose of the tubes after use according to the regulations. Observe the laboratory regulations.

Materials and chemical resistance

The containers are manufactured from polypropylene (PP), polyetheretherketone (PEEK), a thermoplastic elastomer (TPE), and polyethersulfone (PES), respectively polyetheretherketone (PEEK) with PTFE.

The balls are manufactured from either stainless steel AISI 304 or soda-lime glass.

The used plastic materials conform with the FDA (Food and Drug Administration Approval) - rules.

The parts have a good stability against weak acids, chlorides, hypochlorides and many other chemicals.

Technical data

Designation	For disperser unit	Working range	Stator ø	Rotor ø	Circumferential speed at 6,000 min ⁻¹
ST-20...	Tube Drive	2-15 ml	-	-	-
DT-20...	Tube Drive	5-15 ml	26 mm	20 mm	6.5 m/s
BMT-20 ...	Tube Drive	2-15 ml	-	-	-
ST-50...	Tube Drive	15-50 ml	-	-	-
DT-50...	Tube Drive	15-50 ml	26 mm	20 mm	6,5 m/s
BMT-50...	Tube Drive	15-50 ml	-	-	-

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