




AC- 6501xx	<b>CUSTOMER DRAWING</b>	Page 1 of 1
Revision 01	PS Microplate, 96 Well, U-bottom	
	Item-No. 6501xx	



NOT TO SCALE



Cat.- No.	Description	A	B	C	D	E	F	G	H	I
650101	Microplate	127.76 ± 0.2	85.48 ± 0.2	14.2 ± 0.2	∅ 6.94	2.5	10.3	14.29	11.18	9.02
650160	Microplate	127.76 ± 0.2	85.48 ± 0.2	14.2 ± 0.2	∅ 6.94	2.5	10.3	14.29	11.18	9.02
650161	Microplate	127.76 ± 0.2	85.48 ± 0.2	14.2 ± 0.2	∅ 6.94	2.5	10.3	14.29	11.18	9.02
650180	Microplate	127.76 ± 0.2	85.48 ± 0.2	14.2 ± 0.2	∅ 6.94	2.5	10.3	14.29	11.18	9.02
650185	Microplate	127.76 ± 0.2	85.48 ± 0.2	14.2 ± 0.2	∅ 6.94	2.5	10.3	14.29	11.18	9.02

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date April 27, 2005	Date April 27, 2005	Date April 27, 2005	
Date April 27, 2005	Name F. Falter	Name F. Sonntag	Name F. Sonntag	