




AC- 78107x/ 78108x	CUSTOMER DRAWING	Page 1 of 2
Revision 03	384 Well Microplate, PS	 greiner bio-one
	Item-No. 781 073 / 781 074 / 781 075 / 781 076 / 781 077 / 781 079 / 781 080 / 781 086	




NOT TO SCALE



Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date January 22, 2008	Date January 22, 2008	Date January 23, 2008	
Date October 12, 2007	Name F. Falter	Name F. Sonntag	Name F. Sonntag	

AC- 78107x/ 78108x	CUSTOMER DRAWING	Page 2 of 2
Revision 03	384 Well Microplate, PS	
	Item-No. 781 073 / 781 074 / 781 075 / 781 076 / 781 077 / 781 079 / 781 080 / 781 086	

Cat.- No.	Description	A	B	C	D	D'	E	F	G	H	I
781 073	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 074	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 075	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 076	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 077	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 079	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 080	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 086	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date January 22, 2008	Date January 22, 2008	Date January 23, 2008	
Date October 12, 2007	Name F. Falter	Name F. Sonntag	Name F. Sonntag	