



AC- 78109x	<b>CUSTOMER DRAWING</b>	Page 1 of 2
Revision 03	384 Well Cell Culture Microplate, PS, $\mu$ Clear <sup>®</sup>	 <b>greiner bio-one</b>
	Item-No. 781 09x	



NOT TO SCALE



Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date January 22, 2009	Date January 22, 2009	Date January 23, 2009	
Date October 16, 2007	Name F. Falter	Name F. Sonntag	Name F. Sonntag	

AC- 78109x	<b>CUSTOMER DRAWING</b>	Page 2 of 2
Revision 03	384 Well Cell Culture Microplate, PS, µClear®	 greiner bio-one
	Item-No. 781 09x	

Cat.- No.	Description	A	B	C	D	D`	E	F	G	H	I
781 091	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 092	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 093	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 094	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 095	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 096	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 097	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 098	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date January 22, 2009	Date January 22, 2009	Date January 23, 2009	
Date October 16, 2007	Name F. Falter	Name F. Sonntag	Name F. Sonntag	