




AC- 781061/ 7811xx	CUSTOMER DRAWING	Page 1 of 2
Revision 03	384 Well Microplate, PS, Solid Bottom	 greiner bio-one
	Item-No. 781 061/ 781 1xx	




NOT TO SCALE



Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date January 22, 2009	Date January 22, 2009	Date January 23, 2009	
Date October 12, 2007	Name F. Falter	Name F. Sonntag	Name F. Sonntag	

AC- 781061/ 7811xx	CUSTOMER DRAWING	Page 2 of 2
Revision 03	384 Well Microplate, PS, Solid Bottom	 greiner bio-one
	Item-No. 781 061/ 781 1xx	

Cat.- No.	Description	A	B	C	D	D'	E	F	G	H	I
781 061	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 101	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 162	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 165	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 182	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 185	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5
781 186	Microplate	127.76 ±0.15	85.48 ±0.15	14.4 ±0.1	3.7	3.3	2.5	11.5	12.13	8.99	4.5

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date January 22, 2009	Date January 22, 2009	Date January 23, 2009	
Date October 12, 2007	Name F. Falter	Name F. Sonntag	Name F. Sonntag	