




AC- 78208x/ 782180	<b>CUSTOMER DRAWING</b>	Page 1 of 1
Revision 01	1536 Well Cell Culture Microplate CELLSTAR®, PS, HiBase	 greiner bio-one
	Item-No. 782 08x/ 782 180	



Cat.- No.	Description	A	B	C	D	E	F	G	H	I	J
782 080	Microplate	127,76 ±0,1	85,48 ±0,1	10,4 ±0,1	□1,7	□1,53	5	11.01	7.87	2.25	2
782 086	Microplate	127,76 ±0,1	85,48 ±0,1	10,4 ±0,1	□1,7	□1,53	5	11.01	7.87	2.25	2
782 180	Microplate	127,76 ±0,1	85,48 ±0,1	10,4 ±0,1	□1,7	□1,53	5	11.01	7.87	2.25	2

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date April 5, 2006	Date April 6, 2006	Date April 6, 2006	
Date	Name F. Falter	Name F. Sonntag	Name F. Sonntag	