



# Product Specification

**Catalog Number:** BDH9272

**Catalog Description:** BDH POTASS SODIUM TARTRAT

**Grade:** ACS

**Date Current Specification Set:** NOV. 2012

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

Characteristics	Specifications
Ammonium	<= 0.002 %
Assay	99.0 - 102.0 %
Calcium	<= 0.005 %
Chloride	<= 0.001 %
Heavy Metals	<= 5 ppm
Insoluble Matter	<= 0.005 %
Iron	<= 0.001 %
pH (5%, Water) @25.0C	6.0 - 8.5
Phosphate	<= 0.002 %
Sulfate	<= 0.005 %

Signed on behalf of VWR:

Robin Crocker  
Quality Control Manager