



Product Specification

Catalog Number: BDH9276

Catalog Description: BDH SILVER NITRATE ACS

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics	Specifications
Chloride	$\leq 0.0005 \%$
Clarity of Solution	PASS
Copper	$\leq 0.0002 \%$
Free Acid Test	PASS
Iron	$\leq 0.0002 \%$
Lead	$\leq 0.001 \%$
Purity	$\geq 99.0 \%$
Substances not Precipitated by HCl	$\leq 0.01 \%$
Sulfate	$\leq 0.002 \%$

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager