



# Product Specification

**Catalog Number:** BDH9304

**Catalog Description:** BDH STANNOUS CHLORIDE ACS

**Grade:** ACS

**Date Current Specification Set:** NOV. 2012

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

Characteristics	Specifications
<b>TEST</b>	<b>SPECIFICATION</b>
Assay	98.0 - 103.0 %
Calcium	<= 0.005 %
Iron	<= 0.003 %
Lead	<= 0.01 %
Potassium	<= 0.005 %
Sodium	<= 0.01 %
Solubility in hydrochloric acid	PASS
Sulfate (P/F)	PASS

Signed on behalf of VWR:

Robin Crocker  
Quality Control Manager