



# Product Specification

**Catalog Number:** BDH9310

**Catalog Description:** BDH TRICHLOROACETIC ACID ACS

**Grade:** ACS

**Date Current Specification Set:** NOV. 2012

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

Characteristics	Specifications
-----------------	----------------

Chloride	<= 0.002 %
Clarity of Solution	PASS
Heavy Metals	<= 0.002 %
Insoluble Matter	<= 0.01 %
Iron	<= 0.001 %
Nitrate	<= 0.002 %
Phosphate	<= 5 ppm
Purity	>= 99.0 %
Residue after Ignition	<= 0.03 %
Substances Darkened by H2SO4	PASS
Sulfate	<= 0.02 %

Signed on behalf of VWR:

Robin Crocker  
Quality Control Manager