


COMPLETE CONSISTS OF 1 EACH:  
 BOTTLE  
 ORANGE CAP, PP  
 POUR RING, PP



ALL DIMENSIONS ARE NOMINAL  
 DIMENSIONS IN MILLIMETERS UNLESS OTHERWISE NOTED

CORNING INCORPORATED  
 LIFE SCIENCES DIVISION  
 BIG FLATS, NY

ITEM: 1395IS-5L, MEDIA BOTTLE, 5000 mL

ENAMEL COLOR: WHITE

MATERIAL: BOROSILICATE, 33 EXPANSION GLASS

SCALE: 1/2

DATE: 3 OCT 2012

**CORNING**