

Corning® glutagro™ Medium and Supplements

The Best Research Requires The Healthiest Cells.



Corning glutagro medium and supplement combine the quality you expect from Corning with the reagent stability you need. Corning glutagro medium and supplements are the perfect replacement for glutamine.


Cells grown in traditional media supplemented with L-glutamine are susceptible to cytotoxicity. With Corning glutagro medium and supplement, cells receive this needed supplement in the healthiest environment possible.

This modified form of L-glutamine does not degrade into ammonia during storage and incubation, minimizing cytotoxicity and improving culture quality.

Ordering Information

VWR Cat. No.	Corning Cat. No.	Description	Unit Size	Qty/Pk
45001-086	25-015-CI	Corning glutagro medium and supplement	100 mL	1
71002-872	10-101-CV	DMEM with glutagro supplement	500 mL	6
71002-874	10-102-CV	DMEM without sodium pyruvate	500 mL	6
71002-876	10-103-CV	DMEM/F-12 with glutagro supplement	500 mL	6
71002-878	10-104-CV	RPMI 1640 with glutagro supplement	500 mL	6

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.



At Corning, we continuously strive towards improving efficiencies and developing new products and technologies for life science researchers. We have scientists working in Corning R&D labs across the globe, doing what you do every day. Our technical experts understand your challenges and your increased need for high-quality products.

It is this expertise, plus a 160-year legacy of Corning innovation and manufacturing excellence, that puts us in a unique position to be able to offer a beginning-to-end portfolio of high-quality, reliable life sciences consumables.

For additional product or technical information, visit www.vwr.com/corningcellcultureworkflow, call 1.800.932.5000, or contact your VWR representative.

CORNING | **FALCON** | **AXYGEN** | **GOSELIN** | **PYREX**