



## Buy 1 Case of WHEATON Vial Racks Receive 1 FREE Case

### Vial Racks

- Manufactured from polypropylene
- Size of rack depends on vial OD
- Vial OD cannot exceed the well ID dimensions
- Easy to clean and autoclavable
- Alphanumeric indexing indicated on rack
- The racks are sturdy and stackable, even when the vials are in place
- Each well has an opening in the bottom to facilitate drainage

Cat. No.	Description	Well ID (mm)	No. of Wells	Dimensions (L x W x H)	Qty / Case
66022-580	96-Position Rack	9.5	8 deep x 12 wide	18.5 x 12.5 x 2.2cm	5
66022-478	50-Position Rack	12.5	5 deep x 10 wide	19 x 10 x 2.2cm	5
10066-476*	50-Position Rack	12.5*	5 deep x 10 wide	19 x 10 x 2.2cm	5
66022-523	48-Position Rack	15.5	4 deep x 12 wide	26.6 x 9.4 x 2.8cm	5
66022-569	90-Position Rack	17.1	6 deep x 15 wide	33 x 17 x 3cm	5
66022-556	36-Position Rack	23.1	3 deep x 12 wide	32.2 x 9.1 x 2.8cm	5
66022-525	50-Position Rack	28.1	5 deep x 10 wide	33 x 17 x 3cm	5
14007-178	50-Position Rack	30.0	5 deep x 10 wide	33.7 x 17 x 3cm	5

\*This rack is specially designed for use with Freestanding CryoELITE® Cryogenic Vials. The vials lock in the rack for one-handed cap removal.

Purchase your Wheaton Vial Racks by visiting [vwr.com](http://vwr.com) between May 1 and December 31, 2017. To redeem your FREE product: Have your VWR invoice or packing slip ready, visit [vwr.com/promotions](http://vwr.com/promotions), enter the promo code **4985** in the search box and fill in the online form before January 31, 2018. Please allow up to 8 weeks for your free product.

Offer valid in North America. Limited to two redemptions per customer. Free case must be of equal or lesser value. Not for resale. Offer valid while supplies last. All purchases must be on one purchase order.



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.

0517 Lit. No. 230013W Promo Code 4985