

Certificate of Compliance

Customer	VWR International
Customer Catalog #	490001-214
Purchase Order#	45082998072
Product Description	AlumaSeal II™, Custom-Size, Non-Sterile, 100/Box
Excel Catalog #	AF-GBNL-100
Lot#	YF106F
Material:	
Top Coat	N/A
Film	1.4 mil Aluminum Foil
Adhesive	1.1 mil Acrylic
Release Liner	4.6 mil Brown Poly-Coated Paper
Product Requirements	Product is Certified to Meet the Specifications, Dimensions and Tolerances Detailed in Product Drawing# MFG-DWG-112. Product is Certified to Meet the QC Visual & Performance Specifications Detailed in Excel Scientific's QC Procedures. Product is certified free of detectable levels of RNase, DNase and Human/Mouse DNA contamination (PCR Certified).
Packaging	Product has been Packaged in Compliance with BioExpress's Guidelines.
Sterility (if applicable)	Product has been subjected to a gamma radiation level appropriate to achieve 10^{-3} SAL.
Product Shelf-Life:	
Storage Requirements	3 Years if Maintained at 70°F +/-10°; 50% RH +/- 10%
Best if Used By Date	2021-04
Manufacturing Location	USA
Quantity Shipped	500 Boxes
Date of Manufacture	2018-04-16

The Plate-Sealing and Reagent-Handling ExpertsSM

18350 George Blvd. • Victorville, CA 92394 • Tel: 760.246.4545 • Fax: 760.246.5588 • www.excelscientific.com