

Date Prepared: 03/03/2011

Reviewed On: 03/03/2011

### 1 Identification of the substance/mixture and of the company/undertaking

- **Product Identifier**
- **Product Name:** Marine Agar 2216
- **Catalog Number:** 212185
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
 BD Diagnostic Systems  
 7 Loveton Circle  
 Sparks, MD 21152  
 Telephone.: (410) 771 - 0100 or (800) 638 – 8663
- **Information Department:** Technical Service
- **Emergency telephone number:**  
 In case of a chemical emergency, spill, fire, exposure, or accident contact BD Diagnostic Systems (410) 771-0100 or (800)-638-8663, or ChemTrec at (800) 424-9300.


### 2 Composition/information on ingredients

- **Chemical characterization:** Mixture
- **Description:** Mixture consisting of the following components.

- **Dangerous Components:**

7647-14-5	sodium chloride	35.4%
10043-52-4	calcium chloride	3.3%

### 3 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**  
 Xi; Irritant  
 R36/37/38: Irritating to eyes, respiratory system and skin.
- **Classification system:**  
 The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.

---

- **Label elements**
- **Labelling according to EU guidelines:**  
 Observe the general safety regulations when handling chemicals  
 The product has been classified and marked in accordance with regulations on hazardous materials.
- **Code letter and hazard designation of product:**  
 Xi Irritant
- **Hazard-determining components of labelling:**  
 sodium chloride
- **Risk phrases:**  
 36/37/38 Irritating to eyes, respiratory system and skin.
- **Safety phrases:**  
 7 Keep container tightly closed.  
 22 Do not breathe dust.  
 26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.  
 36 Wear suitable protective clothing.

(Contd. on page 2)

Date Prepared: 03/03/2011

Reviewed On: 03/03/2011

**Product Name: Marine Agar 2216**

(Contd. of page 1)

· **NFPA ratings (scale 0-4)**· **HMIS ratings (scale 0-4)**

HEALTH	1	Health = 1
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

**4 First aid measures**

- **General information** No special measures required.
- **After inhalation** Supply fresh air; consult doctor in case of complaints.
- **After skin contact** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact** Rinse opened eye for 15 minutes under running water. Then consult a doctor.
- **After swallowing** If symptoms persist consult doctor.

**5 Firefighting measures**

- **Suitable extinguishing agents**  
CO<sub>2</sub>, ABC multipurpose dry chemical or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Protective equipment:** No special measures required.

**6 Accidental release measures**

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Wipe up with damp sponge or mop.
- **Methods and material for containment and cleaning up:** No special measures required.
- **Reference to other sections** No dangerous substances are released.

**7 Handling and storage**

- **Handling**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Storage**
- **Requirements to be met by storerooms and receptacles:** < 30 °C
- **Information about storage in one common storage facility:** Store away from oxidizing agents.
- **Further information about storage conditions:** None.

**8 Exposure controls/personal protection**

- **Additional information about design of technical systems:** No further data; see Section 7.

(Contd. on page 3)

Date Prepared: 03/03/2011

Reviewed On: 03/03/2011

**Product Name: Marine Agar 2216**

(Contd. of page 2)

- **Components with limit values that require monitoring at the workplace:**  
 The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Personal Protective Equipment**
- **General protective and hygienic measures**  
 The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:**  
 In case of brief exposure, use a chemical fume hood or a NIOSH/MSHA-approved respirator.
- **Protection of hands:**



Chemical resistant gloves (i.e. nitrile, or equivalent).

- **Eye protection:** Safety glasses
- **Body protection:** Protective work clothing (lab coat).

## 9 Physical and chemical properties

- **General Information**
- **Appearance:**
  - Form:** Solid.  
Powder
  - Color:** Beige
- **Odor:** Characteristic
- **Change in condition** Undetermined
- **Melting point/Melting range:** Not determined
- **Boiling point/Boiling range:** Not determined
- **Flash point:** Not applicable
- **Danger of explosion:** Product does not present an explosion hazard.
- **Density:** Not determined
- **Solubility in / Miscibility with**
  - Water:** Soluble

## 10 Stability and reactivity

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Incompatible materials:** Incompatible material: strong oxidizers.
- **Hazardous decomposition products:** Chlorine (Cl)

## 11 Toxicological information

- **Acute toxicity:**
- **Primary irritant effect:**
  - **on the skin:** Irritant to skin and mucous membranes.
  - **on the eye:** Irritating effect.

(Contd. on page 4)

Date Prepared: 03/03/2011

Reviewed On: 03/03/2011

**Product Name: Marine Agar 2216**

(Contd. of page 3)

- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**  
The product shows the following dangers according to internally approved calculation methods for preparations:  
Irritant

**12 Ecological information**

- **Aquatic toxicity:** No further relevant information available.
- **Ecotoxicological effects:**
- **Other information:**  
The ecological effects have not been thoroughly investigated, but currently none have been identified.
- **Additional ecological information:**
- **General notes:** Water danger class 3 (Self-assessment): extremely hazardous for water.

**13 Disposal considerations**

- **Waste treatment methods**
- **Recommendation**  
Smaller quantities can be disposed of with solid waste.  
Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.  
This product is not considered a RCRA hazardous waste.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to state and federal regulations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

**14 Transport information**

- **DOT regulations:**
- **Hazard class:** -
- **Land transport ADR/RID (cross-border)**
- **ADR/RID class:** -
- **Maritime transport IMDG:**
- **IMDG Class:** -
- **Marine pollutant:** No
- **Air transport ICAO-TI and IATA-DGR:**
- **ICAO/IATA Class:** -

**15 Regulatory information**

- **SARA Section 355 (extremely hazardous substances)**  
None of the ingredients is listed.
- **SARA Section 313 (specific toxic chemical listings)**  
None of the ingredients is listed.

(Contd. on page 5)

Date Prepared: 03/03/2011

Reviewed On: 03/03/2011

**Product Name: Marine Agar 2216**

(Contd. of page 4)

· <b>TSCA (Toxic Substances Control Act)</b>	
7647-14-5	sodium chloride
9002-18-0	agar
7786-30-3	magnesium chloride
	sodium sulfate
10043-52-4	calcium chloride
8013-01-2	yeast extract
7447-40-7	potassium chloride
144-55-8	sodium bicarbonate
2338-05-8	Ferric Citrate
7758-02-3	potassium bromide

· **California Proposition 65 - Chemicals known to cause cancer**  
 None of the ingredients is listed.

· **California Proposition 65 - Chemicals known to cause reproductive toxicity for females:**  
 None of the ingredients is listed.

· **California Proposition 65 - Chemicals known to cause reproductive toxicity for males:**  
 None of the ingredients is listed.

· **California Proposition 65 - Chemicals known to cause developmental toxicity:**  
 None of the ingredients is listed.

· **Carcinogenic categories**  
 · **NTP (National Toxicology Program)**  
 None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**  
 None of the ingredients is listed.

· **National regulations**  
 · **Water hazard class:** Water danger class 3 (Self-assessment): extremely hazardous for water.

### 16 Other information

To the best of our knowledge, the information contained herein is accurate. However, neither Becton, Dickinson and Company or any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we can not guarantee that these are the only hazards that exist.

· **Department issuing MSDS:**  
 Environmental, Health & Safety  
 Created by Michael J. Spinazzola  
 · **Contact:** Technical Service Representative