

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

## 1 Identification of the substance/mixture and of the company/undertaking

- **Product Identifier**
- **Product Name:** Tryptic Soy Agar
- **Catalog Number:** 236950
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
 BD Diagnostic Systems  
 7 Loveton Circle  
 Sparks, MD 21152  
 Telephone.: (410) 771 - 0100 or (800) 638 – 8663
- **Information Department:** Technical Service
- **Emergency telephone number:**  
 In case of a chemical emergency, spill, fire, exposure, or accident contact BD Diagnostic Systems (410) 771-0100 or (800)-638-8663, or ChemTrec at (800) 424-9300.


## 2 Composition/information on ingredients

- **Description:** Mixture consisting of the following components.

- **Dangerous Components:**

7647-14-5	sodium chloride	12.22%
-----------	-----------------	--------

## 3 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**  
 Xi; Irritant  
*R36/37/38: Irritating to eyes, respiratory system and skin.*
- **Classification system:**  
*The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.*

---

- **Label elements**

- **Labelling according to EU guidelines:**  
*The product has been classified and marked in accordance with regulations on hazardous materials.*

- **Code letter and hazard designation of product:**



- **Hazard-determining components of labelling:**  
 sodium chloride

(Contd. on page 2)

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: Tryptic Soy Agar**

(Contd. of page 1)

**· Risk phrases:***36/37/38 Irritating to eyes, respiratory system and skin.***· Safety phrases:***7 Keep container tightly closed.**22 Do not breathe dust.**26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.**36 Wear suitable protective clothing.***· NFPA ratings (scale 0-4)**

Health = 1

Fire = 0

Reactivity = 0

**· HMIS ratings (scale 0-4)**

Health = 1

Fire = 0

Reactivity = 0

**4 First aid measures**

- **After inhalation** Seek medical treatment in case of complaints.
- **After skin contact** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact**  
*Rinse opened eye for 15 minutes under running water. Then consult a doctor.*
- **After swallowing** If symptoms persist consult doctor.
- **Information for doctor** Show this product label or this MSDS.

**5 Firefighting measures**

- **Suitable extinguishing agents**  
*CO2, ABC multipurpose dry chemical or water spray. Fight larger fires with water spray or alcohol resistant foam.*
- **Protective equipment:** No special measures required.

**6 Accidental release measures**

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Wipe up with damp sponge or mop.
- **Methods and material for containment and cleaning up:** No special measures required.
- **Reference to other sections**  
*See Section 7 for information on safe handling.*  
*See Section 8 for information on personal protection equipment.*

(Contd. on page 3)

USA

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: Tryptic Soy Agar**

(Contd. of page 2)

See Section 13 for disposal information.

## 7 Handling and storage

- **Handling**
- **Precautions for safe handling** Ensure good ventilation/exhaustion at the workplace.
- **Information about protection against explosions and fires:**  
No special measures required.
- **Storage**
- **Requirements to be met by storerooms and receptacles:** 2 - 30 °C
- **Information about storage in one common storage facility:**  
Store away from oxidizing agents.
- **Further information about storage conditions:**  
Store in cool, dry conditions in well sealed containers.

## 8 Exposure controls/personal protection

- **Additional information about design of technical systems:**  
No further data; see Section 7.
- **Components with limit values that require monitoring at the workplace:**  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Personal Protective Equipment**
- **General protective and hygienic measures**  
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:**  
In case of brief exposure, use a chemical fume hood or a NIOSH/MSHA-approved respirator.
- **Protection of hands:**



Chemical resistant gloves (i.e. nitrile, or equivalent).

- **Eye protection:** Safety glasses
- **Body protection:** Protective work clothing (lab coat).

## 9 Physical and chemical properties

- **General Information**
- **Appearance:**

<b>Form:</b>	Solid.
<b>Color:</b>	Beige
<b>Odor:</b>	Characteristic

(Contd. on page 4)

**Material Safety Data Sheet**  
acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: Tryptic Soy Agar**

(Contd. of page 3)

· <b>pH-value:</b>	7.3+/-0.2
· <b>Change in condition</b>	Undetermined
· <b>Melting point/Melting range:</b>	Not determined
· <b>Boiling point/Boiling range:</b>	Not determined
· <b>Flash point:</b>	Not applicable
· <b>Auto igniting:</b>	Product is not self igniting.
· <b>Danger of explosion:</b>	Product does not present an explosion hazard.
· <b>Density:</b>	Not determined
· <b>Solubility in / Miscibility with Water:</b>	Soluble

**10 Stability and reactivity**

- **Thermal decomposition / conditions to be avoided:**  
No decomposition if used according to specifications.
- **Incompatible materials:** Incompatible material: strong oxidizers.
- **Hazardous decomposition products:** Chlorine (Cl)

**11 Toxicological information**

- **Acute toxicity:**
- **Primary irritant effect:**
  - **on the skin:** Irritant to skin and mucous membranes.
  - **on the eye:** Irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**  
The product shows the following dangers according to internally approved calculation methods for preparations:  
Irritant

**12 Ecological information**

- **Acquatic toxicity:** No further relevant information available.
- **Ecotoxicological effects:**
- **Other information:**  
The ecological effects have not been thoroughly investigated, but currently none have been identified.

(Contd. on page 5)

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: Tryptic Soy Agar**

(Contd. of page 4)

- **Additional ecological information:**
- **General notes:** Not known to be hazardous to water.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**  
 Smaller quantities can be disposed of with solid waste.  
 Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.  
 This product is not considered a RCRA hazardous waste.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to state and federal regulations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

### 14 Transport information

- **DOT regulations:**
- **Hazard class:** -
- **Land transport ADR/RID (cross-border)**
- **ADR/RID class:** -
- **Maritime transport IMDG:**
- **IMDG Class:** -
- **Marine pollutant:** No
- **Air transport ICAO-TI and IATA-DGR:**
- **ICAO/IATA Class:** -

### 15 Regulatory information

- **SARA Section 355 (extremely hazardous substances)**

None of the ingredients is listed.

- **SARA Section 313 (specific toxic chemical listings)**

None of the ingredients is listed.

- **TSCA (Toxic Substances Control Act)**

9002-18-0	agar
9000-71-9	casein
7647-14-5	sodium chloride
8013-01-2	yeast extract
68-04-2	Sodium citrate
127-17-3	pyruvic acid

(Contd. on page 6)

**Material Safety Data Sheet**  
 acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: Tryptic Soy Agar**

(Contd. of page 5)

497-19-8	sodium carbonate
7487-88-9	magnesium sulfate
10043-52-4	calcium chloride
7631-90-5	sodium bisulfite
7646-85-7	zinc chloride
65-22-5	pyridoxal HCL

**· California Proposition 65 - Chemicals known to cause cancer**

None of the ingredients is listed.

**· California Proposition 65 - Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

**· California Proposition 65 - Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

**· California Proposition 65 - Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

**· Carcinogenic categories**
**· NTP (National Toxicology Program)**

None of the ingredients is listed.

**· TLV (Threshold Limit Value established by ACGIH)**

7631-90-5	sodium bisulfite	A4
-----------	------------------	----

**· Product related hazard information:**

The product has been classified and marked in accordance with regulations on hazardous materials.

**· Hazard symbols:**


Xi Irritant

**· Hazard-determining components of labeling:**

sodium chloride

**· Risk phrases:**

36/37/38 Irritating to eyes, respiratory system and skin.

**· Safety phrases:**

7 Keep container tightly closed.

22 Do not breathe dust.

26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

(Contd. on page 7)

USA



**Material Safety Data Sheet**  
acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: Tryptic Soy Agar**

36 Wear suitable protective clothing.

(Contd. of page 6)

### **16 Other information**

*To the best of our knowledge, the information contained herein is accurate. However, neither Becton, Dickinson and Company or any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we can not guarantee that these are the only hazards that exist.*

- **Department issuing MSDS:**  
Environmental, Health & Safety  
Created by Michael J. Spinazzola
- **Contact:** Technical Service Representative

USA